

DIGITAL CENTRIFUGE

DSC-200SD



DSC-200SD

Features:

- * Metal housing is sturdy and durable.
- * Auto-balancing system can extend motor life and reduce noise.
- * With safety interlock, it can prevent the lid from being opened during operation.
- * Forced fan air convection prevents the motor from overheating.
- * Precision digital timer can be set from 0 to 30 minutes with 30 sec. increment.
- * DSC-200SD can offer short runs by pressing the function key.
- * LED digital display enables users to monitor the present speed value or the remaining time.
- * Automatic brake system can offer rapid deceleration.
- * DSC-200SD features a see-through window.
- * The lid will be opened automatically once the rotor stops spinning.
- * Quality certification: ISO 9001, ISO 13485

Specifications for DSC-200SD:

Model	Motor	Timer (Min.)	Display	Weight (kg) (N.W./G.W.)	Overall dimension (mm) (WxDxH)
DSC-200SD	DC motor	0~30 min. digital timer	LED digital display for time and speed	14.5/16.0	320x420x260

Power supply: AC 110V 60Hz (220V 50/60Hz available)

Specifications of rotor for DSC-200SD:

Model	Speed (rpm)	Type	Max. RCF (g)	Capacity	Standard accessory
UAR-1508	300~6,000	45 degree angle	4,145	8 x (5~15)ml	(5~7)ml and (10~15)ml tube adapters

Note: Tube adapter for (1.5~2.0)ml is optional accessory

Easy order information:

Order No.	Item
DSC-N1508SD	Centrifuge DSC-200SD with rotor UAR-1508
CTA15	Tube adapter for (10~15)ml tube; inner size: Diam. 17.5 x L98mm
CTA7	Tube adapter for (5~7)ml tube; inner size: Diam. 13.5 x L69mm
CTA2	Tube adapter for (1.5~2.0)ml tube; inner size: Diam. 11.3 x L34mm



Rotor UAR-1508

Tube adapters



CTA15

CTA7

CTA2



DIGISYSTEM LABORATORY INSTRUMENTS INC.

NO. 1, LANE 393, FU DE 1 ROAD, XIZHI DIST., NEW TAIPEI CITY, TAIWAN
TEL: 886-2-26941738 . 26942339 . 26945547 FAX: 886-2-26941695
E-MAIL: digilab@ms23.hinet.net INTERNET: www.digisystemlab.com
We reserve the rights to change the specifications without notice.

V.201301